

GEL ROCKER

Manufacturer of
Life Science Research
Equipment




GeM
Government
e Marketplace



EPS offers advance technology Gel Rocker are known for their durability and Rockers are commonly used for staining and de-staining gels after electrophoresis, hybridization, washing, blotting, Cell Culture and gentle mixing other chemical mixing procedures.(See-Saw Movement)

Features

- Variable oscillation model with maximum 60 oscillations per minute.
- Fine adjusting facility for setting the parameters
- MOC: GI Steel treated and powder coated to prevent corrosion.
- Maximum capacity can withstand upto 10 kgs.
- Auto restart set program in case of power failure.
- Timer - 1 Min to 99 hr. 59 Min & Continues Mode
- User Guide



Technical Specifications & Ordering Information: Single Type

MODEL	EPS/GR MONO	EPS/GR MONO LCD
Top platform in mm	250 X 250	250 X 250
Speed	0 to 60 Oscillation	0 to 60 Oscillation
Timer	Nil	1 Min to 99hr .59Min & Continues Mode
Control	Knob	Soft Membrane key
Display	-	LCD (Timer and RPM)

Dual Type

MODEL	EPS/GR TWIN	EPS/GR TWIN LCD
Top platform in mm	2 x (250 X 250)	2 x (250 X 250)
Speed	0 to 60 Oscillation	0 to 60 Oscillation
Timer	Nil	1 Min to 99hr .59Min & Continues Mode
Control	Knob	Soft Membrane key
Display	-	LCD(Timer and RPM)



Madras University



The Leader in Immunotherapy



SRI RAMACHANDRA
UNIVERSITY
Chennai



Our Dealers Network

Educational Institutions ◀

Pharma Industries ◀


Research Organizations ◀

Food & Drug Testing Laboratories ◀

Diagnostic Labs ◀

Aquaculture Labs ◀

Environmental & Water Treatment Units ◀



Authorized Dealer / Distributor

EPS Biosolutions

Registered Office:

Plot 4B, Annai Abode Apartment, 1st Cross Street,
Chitti Babu Nagar, Pallikaranai, Chennai - 600 100.
Tamilnadu. India. ☎ : 044-4207 2939.

✉ : info@epsbiosolutions.com | sales@epsbiosolutions.com

🌐 : www.epsbiosolutions.com

Factory :

Plot No: 201-B, 9th Main Road,
Ram Nagar, South Extension,
Pallikaranai, Chennai-600100, Tamil Nadu, India

☎ : 044-42072939 ✉ : info@epsbiosolutions.com

🌐 : www.epsbiosolutions.com