

UV TRANSILLUMINATOR

Manufacturer of
Life Science Research
Equipment



GeM
Government
e Marketplace



LUMINA

LUMINA high quality UV Transilluminators are designed for viewing DNA/RNA in agarose gels stained with ethidium bromide EB or other dyes.

Salient Features/USP

- Available with different range of wave lengths 254 nm, 312 nm and 360 nm
- The high intensity imported UV tubes are used for high illumination
- Imported UVG long-life filters are rated for 5000 hours of use
- The imported UV filter for high end models
- High/low intensity soft membrane key with LCD Display
- The cent percent UV resistant Polycarbonate shield for safety purpose
- The Nylon hinges provided with shield and open it as long as it can



UV Model

Technical Specifications & Ordering Information

MODEL	EPS/UV15365	EPS/UV15312	EPS/UV20365	EPS/UV20312	EPS/UV365A	EPS/UV15	EPS/UV20
Wavelength	365nm	312nm	365nm	312nm	365nm	365/312/254nm	365/312/254nm
Illumination Area(cm)	15x15	15x15	20x20	20x20	20x15	15x15	20x20
Light source	4x4 Watts	4x8 Watts	6x8 Watts	6x8 Watts	4x8 Watts	6x8 Watts	6x8 Watts
LCD Display	-	-	Optional	Optional	-	-	Optional
Exhaust Fan	Yes						

UV-VIS Model

Technical Specifications & Ordering Information:


MODEL	EPS/UVVIS15365	EPS/UVVIS15312	EPS/UVVIS20365	EPS/UVVIS20312	EPS/UVVIS365A	EPS/UVVIS15D	EPS/UVVIS20D
Wavelength	365nm	312nm	365nm	312nm	365nm	365/312/254nm	365/312/254nm
Illumi. Area (cm)	15x15	15x15	20x20	20x20	20x15	15x15	20x20
Light source	4x4 Watts	4x8 Watts	6x8 Watts	6x8 Watts	4x8 Watts	6x8 Watts	6x8 Watts
LCD Display	-	-	Optional	Optional	-	-	Optional
Exhaust Fan	Yes						



Madras University



SRI RAMACHANDRA
UNIVERSITY
Chennai



Our Dealers Network

Educational Institutions ◀

Pharma Industries ◀


Research Organizations ◀

Food & Drug Testing Laboratories ◀

Diagnostic Labs ◀

Aquaculture Labs ◀

Environmental & Water Treatment Units ◀



Authorized Dealer / Distributor

EPS Biosolutions

Registered Office:

Plot 4B, Annai Abode Apartment, 1st Cross Street,
Chitti Babu Nagar, Pallikaranai, Chennai - 600 100.
Tamilnadu. India. ☎ : 044-4207 2939.

✉ : info@epsbiosolutions.com | sales@epsbiosolutions.com

🌐 : www.epsbiosolutions.com

Factory :

M/s EPS Biosolutions,
Plot No: 201-B, 9th Main Road, Ram Nagar,
South Extension, Pallikaranai,
Chennai - 600100, Tamil Nadu, India
T : 044-42072939 E : info@epsbiosolutions.com
W : www.epsbiosolutions.com