

WATER BATH

Manufacturer of
Life Science Research
Equipment



GeM
Government
e Marketplace



WATER BATH

Water Baths are designed manufactured to comply requirements of industries. Water baths are suitable for various sizes of flasks supplied with concentric rings Holes and racks Serological . Water Baths are fully double wall construction and the gap between the walls insulated with special grade glass wool to eliminate the heat loss to minimum. U Shaped heaters are used.



WATER BATH (HOLES)

Technical Specifications & Ordering Information:

MODEL	EPS/WB-R6	EPS/WB-R9	EPS/WB-R12
Dimension (mm)	300x250x105	300x300x105	405x300x105
No. of Holes Width	6	9	12
Volume Litres	9	12	15
MOC	Outer :- Mild Steel (Powder coated) , Inner :- Stainless steel		
Control	Thermostat/Microprocessor Digital controller (Optional)		
Heaters	'U' Shaped heater		
Power	1.5 kw		
Temperature	Ambient +5°C to 100°C		

WATER BATH (SEROLOGICAL)

EPS Serological Water Baths are fully double wall construction and the gap between the walls insulated with special grade glass wool to eliminate the heat loss to minimum. 'U' Shaped heaters are used.



Technical Specifications & Ordering Information:

MODEL	EPS/WB-S2	EPS/WB-S4	EPS/WB-S6	EPS/WB-S8
Dimension (mm)	300x250x175	350x300x175	455x300x175	605x300x175
Volume Width	15	20	25	32
MOC	Outer :- Mild Steel (Powder coated) , Inner :- Stainless steel			
Control	Thermostat/Microprocessor Digital controller (Optional)			
Heaters	'U' Shaped heater			
Power	1.5 kw			2.0 kw
Temperature	Ambient +5°C to 100°C			



Madras University



The Leader in Immunotherapy



SRI RAMACHANDRA
UNIVERSITY
Chennai



National Institute of Animal Biotechnology



Making Lives Happier



INDIAN COUNCIL OF
MEDICAL RESEARCH
Established in 1952



Our Dealers Network

Educational Institutions ◀

Pharma Industries ◀


Research Organizations ◀

Food & Drug Testing Laboratories ◀

Diagnostic Labs ◀

Aquaculture Labs ◀

Environmental & Water Treatment Units ◀



EPS Biosolutions

Registered Office:

Plot 4B, Annai Abode Apartment, 1st Cross Street,
Chitti Babu Nagar, Pallikaranai, Chennai - 600 100.
Tamilnadu. India. ☎ : 044-4207 2939.

✉ : info@epsbiosolutions.com | sales@epsbiosolutions.com

🌐 : www.epsbiosolutions.com

Factory :

Plot No: 201-B, 9th Main Road,
Ram Nagar, South Extension,
Pallikaranai, Chennai-600100, Tamil Nadu, India
☎ : 044-42072939 ✉ : info@epsbiosolutions.com

🌐 : www.epsbiosolutions.com