

SUBMARINE GEL ELECTROPHORESIS (DUAL / TWIN DUAL)

Manufacturer of
Life Science Research
Equipment



GeM
Government
e Marketplace



SUBMARINE GEL ELECTROPHORESIS (DUAL / TWIN DUAL)


EPS horizontal gel electrophoresis offers many advantages for nucleic acid separation and remains today one of the mainstays of molecular biology. EPS range offers the most comprehensive and versatile range of horizontal gel electrophoresis units currently available for low- and high-throughput DNA and RNA applications.

Features



- Color-coded and height adjustable(optional) combs-easily identify comb thickness at a glance and control well depth
- Can run up to 4 no. of gel running trays in a single time
- Colour background-for easy well detection when loading
- Compact tank-reduces the buffer volume
- Electrodes can be detachable and replaced by many times if needed.
- UV-transparent running tray-allows the user to image the gel without risk of damage due to handling
- Side handles-for safe and easy transportation around the laboratory

Complete Unit Includes

- Buffer tank with safety lid
- Electrodes (Positive and negative)
- Gel running tray, combs
- (Gel casting Assembly-Optional)
- Power cord
- User Guide



Submarine Gel Electrophoresis (Dual)



Submarine Gel Electrophoresis (Twin Dual)

Technical Specifications & Ordering Information :

MODEL	EPS/MINI D SUB 10X7	EPS/MIDI D SUB 13X10	EPS/MAXI D SUB 20X20	EPS/MAXI SUB 26X20 D-TWIN
Buffer Capacity	1.2L	2L	3.5L	4L
Gel size(cm)	10X7	13X10	20X20	26X20
No. of wells	8	15 & 18	26 & 30	26 & 51
Max samples	32	72	120	616
Electrode Diameter (mm)	0.3	0.3	0.3	0.3



Our Dealers Network

Educational Institutions ◀

Pharma Industries ◀


Research Organizations ◀

Food & Drug Testing Laboratories ◀

Diagnostic Labs ◀

Aquaculture Labs ◀

Environmental & Water Treatment Units ◀



EPS Biosolutions

Registered Office:

Plot 4B, Annai Abode Apartment, 1st Cross Street,
Chitti Babu Nagar, Pallikaranai, Chennai - 600 100.
Tamilnadu. India. ☎ : 044-4207 2939.

✉ : info@epsbiosolutions.com | sales@epsbiosolutions.com

🌐 : www.epsbiosolutions.com

Factory :

No. 28/17, Bajanai Kovil Street,
Off - Orandiamman Kovil Street,
Velachery, Chennai - 600 042. Tamil Nadu, India

☎ : 044-42072939 ✉ : info@epsbiosolutions.com

🌐 : www.epsbiosolutions.com