



# MUFFLE FURNACE



Manufacturer of  
Life Science Research  
Equipments



EPS Muffle Furnace offer temperature range starts 900°C, 1100°C, 1200°C, 1300°C & up to 1500°C and temperature control options to meet your application needs. Embedded or open heating elements designed to keep samples safe while maintaining reliable temperature uniformity.

Furnaces in laboratory and industrial settings are used every day for a wide variety of simple and technical applications. Choose from a wide offering to accommodate your applications needs, which may include.

- Metal treatment
- Ashing
- Water treatment
- Pottery Designs with safety in mind

## Temperature:

- Controlled with automatic Digital temperature indicating controller with Cr/A1. Thermocouple

## Construction:

- Outer chamber made up of thick cold rolled mild steel with powder coated for corrosive resistance.
- Chamber is surrounded by ceramic fibre blankets on all the four sides and the rear to reduce heat losses to a minimum.

## Doors:

- Rectangular in shape placed horizontally with opening in the front
- Sturdy handle provides comfortable opening and closing.
- Doors open 180° in order to provide enough access to the heating chamber.



## Technical Specifications & Ordering Information:

MODEL	EPS/MF-100	EPS/MF-125	EPS/MF-150	EPS/MF-200	EPS/MF-300
Inner Chamber (mm)	100x100x250	125x125x250	150x150x300	200x200x300	300x300x300
Load	3.5 KW	4.0 KW	4.5 KW	5.0 KW	6.0 KW
Temperature Range	900°C, 1000°C, 1100°C, 1200°C, 1300°C up to 1500°C				
Temperature Accuracy	+ / - 2°C				
Input Power Supply	230 Volts AC, Single phase				



भारत सरकार Government of India सूक्ष्म, लघु एवं मध्यम उद्यम मंत्रालय Ministry of Micro, Small and Medium Enterprises	
<b>UDYAM REGISTRATION CERTIFICATE</b>	
UDYAM REGISTRATION NUMBER	UDYAM-TN-42-613876
NAME OF ENTERPRISE	MS EPS BIOSOLUTIONS
TYPE OF ENTERPRISE *	MICRO (Based on FY 2019-20)
MAJOR ACTIVITY	MANUFACTURING
SOCIAL CATEGORY OF ENTREPRENEUR	OBC
NAME OF UNIT(S)	1. EPS BIOSOLUTIONS
OFFICIAL ADDRESS OF ENTERPRISE	Flat No. 4, 1st Floor, 1st Cross Street, Chittibabu Nagar, 1st Cross Street, Chittibabu Nagar, Chennai, Chengalpattu, Tamil Nadu, 600100
DATE OF INCORPORATION / REGISTRATION OF ENTERPRISE	16/11/2008
DATE OF COMMENCEMENT OF PRODUCTION/BUSINESS	16/11/2008
NATIONAL INDUSTRY CLASSIFICATION CODE(S)	31.20 - Other manufacturing
DATE OF UDYAM REGISTRATION	06/06/2022

\* In case of graduation (upward movement) of status of an enterprise, the benefits of the Government Schemes will be available as per the provisions of Notification No. S.O. 2186/2018 dated 26.08.2018 issued by the MoSME.

This certificate is a representation of the status of an enterprise as on the date of registration. It is not a guarantee of the status of an enterprise. The status of an enterprise may change over time. The status of an enterprise may be verified by visiting the UDYAM portal.

For any assistance, you may contact:

1. District Industries Centres: CHENNAI (TAMIL NADU)

2. MSME DE: CHENNAI (TAMIL NADU)

Visit : [www.msme.gov.in](http://www.msme.gov.in) | [www.dicmsme.gov.in](http://www.dicmsme.gov.in) | [www.champions.gov.in](http://www.champions.gov.in)

BE A CHAMPION with the Ministry of MSME

भारत सरकार Government of India व्यापार चिह्न रजिस्ट्री Trade Marks Registry व्यापार चिह्न अधिनियम, 1999 Trade Marks Act, 1999	
व्यापार चिह्न नाम / Trade Mark No.	3938903
व्यापार चिह्न प्रकार / Trade Mark Class	9
व्यापार चिह्न विवरण / Trade Mark Description	04th day of April, 2019
व्यापार चिह्न रजिस्ट्री / Trade Marks Registry	MUMBAI

Certified that Trade Mark / a representation is attached herewith, has been registered in the name(s) of - EPS BIOSOLUTIONS, No. 83, F4 VIGNESHWARA APARTMENT, 8th CROSS STREET, RAJALAKSHMINAGAR, VELACHERY, CHENNAI 600 042, K. GOPAL ANTHA, PARTNERSHIP FIRM, (Partnership Firm)

In Class 9 Under No. 3938903 as of the date 06 September 2018 in respect of

Goods and service as annexed

Trade Mark as annexed


व्यापार चिह्न रजिस्ट्री / Trade Marks Registry MUMBAI

व्यापार चिह्न रजिस्ट्री / Trade Marks Registry MUMBAI

## OUR CLIENTS



## Our Dealers Network



**Educational Institutions** ◀

**Pharma Industries** ◀

**Research Organizations** ◀

**Food & Drug Testing Laboratories** ◀

**Diagnostic Labs** ◀

**Aquaculture Labs** ◀

## EPS Biosolutions

### **Registered Office:**

Plot 4B, Annai Abode Apartment, 1st Cross Street,  
Chitti Babu Nagar, Pallikaranai, Chennai - 600 100.

Tamilnadu. India. ☎: 044-4207 2939.

🌐 info@epsbiosolutions.com | sales@epsbiosolutions.com

✉ www.epsbiosolutions.com

### **Factory:**

No. 28/17, Bajanai Kovil Street,  
Velachery, Chennai - 600 042. Tamil Nadu, India

Off - Orandiamman Kovil Street,

☎ 044-42072939 🌐 info@epsbiosolutions.com

✉ www.epsbiosolutions.com